



# DIRECTORY OF SERVICES

## TOXICOLOGY SCREEN

LIS Test #	Test Name	CPT Code	Specimen Type	Device	Transport Temperature	RT	FG	FZ	Preferred Volume	Min. Specimen Volume	Notes
6AM	6-AM (Heroin Metabolite)	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
AMPH	Amphetamines	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
BARB	Barbiturates	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
BENZ	Benzodiazepines	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
BUP	Buprenorphine	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
THC	Cannabinoid	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
COCM	Cocaine	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
XTC	Ecstasy	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
ETOH	Ethyl Alcohol	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
ETG	Ethyl Glucuronide	80302	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
FENT	Fentanyl	80302	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
METH	Methadone	80301	Urine	Urine Cup	RT or FG	7 d	30 d	1 y	2 mL	1 mL	Boric Acid tube cannot be utilized
EDDP	Methadone Metabolite	80301	Urine	Urine Cup	RT or FG	7 d	30 d	2 y	2 mL	1 mL	Boric Acid tube cannot be utilized
OPI	Opiates	80301	Urine	Urine Cup	RT or FG	7 d	30 d	3 y	2 mL	1 mL	Boric Acid tube cannot be utilized
OXY	Oxycodone	80301	Urine	Urine Cup	RT or FG	7 d	30 d	4 y	2 mL	1 mL	Boric Acid tube cannot be utilized
PCP	Phencyclidine	80301	Urine	Urine Cup	RT or FG	7 d	30 d	5 y	2 mL	1 mL	Boric Acid tube cannot be utilized
TCA	Tricyclic	80301	Urine	Urine Cup	RT or FG	7 d	30 d	6 y	2 mL	1 mL	Boric Acid tube cannot be utilized
CREAT-DOA	Validity-Creatinine	82570	Urine	Urine Cup	RT or FG	7 d	30 d	6 y	2 mL	1 mL	Boric Acid tube cannot be utilized
PH	Validity-pH	83986	Urine	Urine Cup	RT or FG	7 d	30 d	6 y	2 mL	1 mL	Boric Acid tube cannot be utilized
GRAV	Validity-Specific Gravity	81003	Urine	Urine Cup	RT or FG	7 d	30 d	6 y	2 mL	1 mL	Boric Acid tube cannot be utilized
OXID	Validity-Oxidant	84311	Urine	Urine Cup	RT or FG	7 d	30 d	6 y	2 mL	1 mL	Boric Acid tube cannot be utilized



# DIRECTORY OF SERVICES

## TOXICOLOGY CONFIRMATION

LIS Test #	Test Name	Compounds Detected	CPT Code	Specimen Type	Device	Transport Temperature	RT	FG	FZ	Preferred Volume	Min. Specimen Volume
6MO	6-AM (Heroin)	6-Acetylmorphine	80356	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
EOHG	Alcohol Biomarker	Ethyl beta-glucuronide	80321 (1-2)	Urine	Urine Cup	RT or FG	7d	7d	7d	2 mL	1 mL
ETS		Ethyl sulfate		Urine	Urine Cup	RT or FG	7d	7d	7d	2 mL	1 mL
MITG	Alkaloids, NOS	Mitragynine	80323	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
KAVA		Kavain (Kava)		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
AMP	Amphetamine	Amphetamine	80325 (3-4)	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
METHA		Methamphetamine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MEPHE		Mephedrone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
AMIT	Antidepressant-TCA	Amitriptyline	80336 (3-5)	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
DESIP		Desipramine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
DOXE		Doxepin		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
IMIP		Imipramine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NAMIT		Nortriptyline		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
BUPROP	Antidepressant-SSRI	Bupropion		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
BUSP		Buspirone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
FLUOX		Fluoxetine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
DESME		O-Desmethylvenlafaxine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
SERT		Sertraline		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
VEN		Venlafaxine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
CLOZ	Antipsychotics	Clozapine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NDESME		N-Desmethylozapine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
BUTB	Barbiturate Class	Butalbital	80345	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
PBARB		Phenobarbital		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
SECO		Secobarbital		Urine	Urine Cup	RT or FG	4d	10d	30d	3 mL	2 mL

7-AM	Benzodiazepine Class	7-Amino-Clonazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
AHY		Alpha hydroxyalprazolam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
AZO		Alprazolam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
DIAZ		Diazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
FLU		Flurazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
LORAZ		Lorazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MIDA		Midazolam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NDIAZ		Nordiazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
OXAZ		Oxazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
TEM		Temazepam		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
BUPR	Buprenorphine	Buprenorphine	80348	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NORBUP		Norbuprenorphine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
THC-C	Cannabinoids; Natural	Delta-9-THC-COOH	80349	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
W250 PEN	Cannabinoids; Synthetic	JWH 250 5OH Pentyl		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
BENZO	Cocaine	Benzoylecgonine	80353	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
FEN	Fentanyl	Fentanyl	80354	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NORFEN T		Norfentanyl		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
CARF		Carfentanil		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NORC		Norcarfentanil		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
GAB	Gabapentin	Gabapentin	80355	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
DLSD	Hallucinogen	D-Lysergic Acid Diethylamide		Urine	Urine Cup	RT or FG	4d	10d	30d	3 mL	2 mL
KETA	Ketamine	Ketamine	80357	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NKETA		Norketamine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MTD	Methadone	Methadone	80358	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
ED		EDDP		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MDA	MDMA	MDA	80359	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MDEA		MDEA		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MDMA		MDMA		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
RA	Methylphenidate	Ritalinic Acid	80360	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL

CODE	Opiate	Codeine	80361 (1 or more)	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
HCOD		Hydrocodone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
HMORPH		Hydromorphone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NHCOD		Norhydrocodone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MORPH		Morphine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NALT	Opioid and Opiate Analog	Naltrexone	80363 (3-4)	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
6BETA		6-beta Naltrexol		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
136D		Dextromethorphan		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
W18		W18		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
NALO		Naloxone		Urine	Urine Cup	RT or FG	4d	10d	30d	3 mL	2 mL
NORME	Normeperdine	Urine	Urine Cup	RT or FG	4d	10d	30d	4 mL	3 mL		
NOXY	Oxycodone	Noroxycodone	80365	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
OXYC		Oxycodone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
OMORPH		Oxymorphone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
PC	Phencylidine (PCP)	Phencyclidine (PCP)	83992	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
ZOL	Sedative Hypnotic	Zolpidem	80368	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
ZOLCA		Zolpidem Carboxylic Acid		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
CYCLO		Cyclobenzaprine		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MEPRA		Meprobamate		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
FLA	Stimulant-Synthetic	alpha-PVP (Flakka)	80371	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
MDPV		MDPV		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
METHY		Methylone		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
TRAM	Tramadol	Tramadol	80373	Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL
O-DES		O-Desmethyltramadol		Urine	Urine Cup	RT or FG	4d	10d	30d	2 mL	1 mL